

## Aspire Lifecomfort Acute Care Mattress



The LIFECOMFORT® ACUTE CARE mattress has three foam layers including a castellated patient surface providing users with immersive, pressure relieving support and air flow for enhanced skin microclimate. The U-shaped base and sidewall structure provides support to assist with patient transfers as well as creates roll resistance at the sleep surface edges. Fitted with superior, hospital-grade polyurethane coated covers designed to withstand the rigors of usage and cleaning requirements in Home Care, Aged Care and up to Acute Care environments.

- **Durable Hospital Grade Foam** that is treated for antimicrobial and antibacterial properties.
- **Breathable, 4-Way Stretch Waterproof Cover** to reduce shear, friction and maintain skin microclimate.
- **Medical Grade Waterfall Cover** is fully-welded over the zipper to maximise infection control.
- **Asset Management Ready** with unique batch and serial numbers printed on all covers.
- **Carry Handle** Integrated handles for safe manual handling.
- **Hospital Grade** - Suitable for bleach concentrations up to 10,000ppm.
- **4-sided Zipper** - Easily removed for effective cleaning.
- **4-way Stretch** - For shear mitigation and friction reduction.
- Immersive surface contours and cradles around the user for personalised pressure relief and comfort.
- Machine Washable Covers: Up to 95O
- **Foam Layers:**
  - IMMERSIVE CASTELLATED SURFACE - Castellated cells offer greater airflow and transitional support
  - COMFORT LAYER - Visco-elastic comfort and envelopment
  - SAFETY & SUPPORT BASE - Strong U-shaped base for maximum durability

and assistive support

- INTEGRATED SIDE WALLS - Strong rigid side layers for transfer efficiency and safety
- **Applications** - Pressure Care: Medium - High
  - Hospital
  - Aged Care
  - Home Care
  - Palliative Care

## Warranty Information

- 5 yrs Foam
- 4 yrs Cover

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
UK-BEM046615	Lifecomfort Acute Care Mattress - Single	6'6"	34"	6"	9 - 39.5 st	18.8 lb
UK-BEM046616	Lifecomfort Acute Care Mattress - King Single	6'8"	3'5"	6"	9 - 39.5 st	18.8 lb