

Aspire Lifecomfort Bed Rail Protector



- Patient Protection - Protects patients from injury by preventing limbs from protruding through bed rails as well as any contact
- Velcro Fastening Straps - Easy to apply and remove
- Functionality - Does not impede routine bed rail functionality
- Easy Clean PU Fabric - Covers are easily wipe cleaned to assist with infection control management
- Sold in pairs and available in 2 sizes for optimal fit
- Washable up to 95°C with water and mild detergent

Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width
UK-BEA010775	Lifecomfort Bed Rail Protector - Small	4'1"	26"
UK-BEA010776	Lifecomfort Bed Rail Protector - Large	5'2"	38"