

Aspire STAR Mover



- The Aspire Star Mover offers a safe and comfortable solution for patients and carers doing sit to stand transfers.
- The compact and versatile design, is perfect for Homecare, Institutional and Aged Care environments.
- Wide seat entry suits a larger range of body shapes.
- Low footplate height makes patient entry and exit easier.
- Large 3.9' rear locking castors have low rolling resistance and offer a smooth ride during transportation.
- Narrow base improves access in tight spaces.
- Height adjustable knee support ensures correct positioning for a larger patient height range.
- Easy access seat pad handles provide reach-around access to carers.
- Soft touch grab bar provides a comfortable patient support for improved stability and confidence during transfers.
- Width adjustable legs with tapered close position allow access between and, or around chair legs.

Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
UK-LSS390580	Aspire STAR Mover	36"	24" – 35"	3'5"	29 st