

Aspire Flex Adjustable Seat Walker



- Super lightweight making it easy to manoeuvre and easy to transport.
- Available in 6' (Lilac & Champagne) and 8' (Blue & Black) wheel variations.
- Robust aluminium frame for improved durability and stability.
- Height adjustable seat height to suit most users.
- Soft padded seat for improved seating comfort.
- Adjustable and lockable brake handles.

Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Width	Seat To Floor Height	Seat To Footrest
UK-WAF709660CH	Aspire Flex Adjustable Seat Walker 6" Wheel - Champagne	28.7"	22.4"	28.7 - 32.6"	21 st	13.8 lbs	14.1"	16.9 - 20.4"	16.9 - 20.4"
UK-WAF709660PU	Aspire Flex Adjustable Seat Walker 6" Wheel - Lilac	28.7"	22.4"	28.7 - 32.6"	21 st	13.8 lbs	14.1"	16.9 - 20.4"	16.9 - 20.4"
UK-WAF709670BK	Aspire Flex Adjustable Seat Walker 8" Wheel - Black	32.2"	23"	30.7 - 38.1"	21 st	16 lbs	14.1"	18.7 - 22.4"	18.7 - 22.4"
UK-WAF709670BL	Aspire Flex Adjustable Seat Walker 8"	32.2"	23"	30.7 -	21 st	16 lbs	14.1"	18.7 - 22.4"	18.7 - 22.4"

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Width	Seat To Floor Height	Seat To Footrest
	Wheel - Blue			38.1"					