

## KeraSit2Sit Accessories



- Kera Site2Sit Accesories
- **Back Strap** - Specially designed Kera sit2sit footboxes, which are easily fitted, to hold the user's feet in place. Sold as a pair, complete with footbox straps and toe guards.
- **Thighstrap** - The Thighstrap is easily added to the Kera sit2sit and reduces the tightness required from the Backstrap.
- **Riser Blocks** - Add these EVA riser blocks (2.36 ' or 4.33 ') when the user's feet don't reach the footplates. Designed to be used with Kera sit2sit footboxes.
- **Heel Straps** - Made for use with Kera Footboxes to secure patients feet during transfer.
- **Foam Packs** - 4 x foam packers designed to fit behind the Kera sit2sit Knee Support Covers.

Code	Display Name
UK-LSS398300	Back Strap - Kera Sit2Sit
UK-LSS398310	Thigh Straps - Kera Sit2Sit
UK-LSS398320	Foot Box - Kera Sit2Sit
UK-LSS398330	Heel Straps - Kera Sit2Sit
UK-LSS398340	Foam Packs - Kera Sit2Sit