

4H415G Tilt in Space Shower Commode



- Robust, lightweight stainless steel design
- Gas assisted tilt in space design
- Polyurethane backrest which is removable for cleaning purposes
- Padded swing away armrests
- Epoxy powder coated frame
- Height adjustable swing-away legrests with calf padding
- 5' tracking swivel castors
- Anti-tippers as standard
- Closed padded seat as standard
- Plastic pan/bowl runners

Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

Code	Display Name	Safe Working Load	Seat Width	Backrest Height
UK-BTC069470	4H415G Tilt in Space Shower Commode - 22" 28 st	22"	-	
UK-BTC069500	4H415G Tilt in Space Shower Commode - 19" 28 st	19"	35"	