

## Aspire Evoke 2 Matte Black Wheelchair



- The Evoke 2 is designed for the user who wants a lightweight, highly adjustable manual wheelchair
- The Contoured Padded Overlay is a stylish recommended accessory for the Aspire Evoke 2. The contoured padded overlay offers the user increased comfort and is ideal for those who are in their wheelchair for longer periods of time. The mild contouring assists in positioning and pressure distribution, while the fabric allows air movement, increasing breathability. It is easily removed to allow washing when required
- Swing-in / out, height adjustable and removable foot rests – this facilitates user transfers and allows optimal positioning for the changing needs of the user
- Strong, Nylon upholstery with a comfortable cushion and padded backrest. For increased comfort, a padded overlay is available as an accessory
- Depth adjustable seating with rear angle adjustability – for more accurate positioning, customisability and increased user comfort
- Quick release 24" solid rear tyres - the wheels are easily removed to reduce the weight of the chair for transport and storage.
- Crash Tested – to AS/NZ3616.196

### Warranty Information

- 2 Years Frame 1 Year Components

Code	Display Name	Total Width	Safe Working Load	Colour	Seat Depth	Seat Width	Backrest Height
UK-MWS449505	Aspire Evoke 2 Matte Black Wheelchair Standard 14" 22"	24 st		Matte Black	12" – 18"	14"	17" – 19"

Code	Display Name	Total Width	Safe Working Load	Colour	Seat Depth	Seat Width	Backrest Height
UK-MWS449515	Aspire Evoke 2 Matte Black Wheelchair Standard 16" 24"	24"	24 st	Matte Black	12" – 18"	16"	17" – 19"
UK-MWS449525	Aspire Evoke 2 Matte Black Wheelchair Standard 18" 26"	26"	24 st	Matte Black	12" – 18"	18"	17" – 19"
UK-MWS449535	Aspire Evoke 2 Matte Black Wheelchair Standard 20" 28"	28"	24 st	Matte Black	16" – 22"	20"	17" – 19"