

## Aspire Shower Chair



- Safe Working Load (SWL) of 31 st tested to AS/NZS39737:2009
- Contoured plastic seat with drainage holes and windows - Provides postural support, improves drainage and reduces skin adhesion
- Height adjustable legs - Easy push button height adjustment for changing needs
- Additional height increments on one leg - Accommodates uneven surfaces
- Raised armrests - Assist the user transferring on and off the chair
- All frames are individually treated for corrosion resistance
- Durable non-slip feet - Provide excellent grip and stability in wet conditions and a stainless steel washer within enhances durability

### Warranty Information

- Corrosion: 5 Years
- Components: 1 Year



| Code         | Display Name                           | Total Width | Safe Working Load | Unit Weight | Seat Depth | Seat Width | Seat To Floor Height | Seat To Footrest |
|--------------|--|-------------|-------------------|-------------|------------|------------|----------------------|------------------|
| UK-BTS118000 | Shower Chair with Arms - Treated Steel | 20"         | 31.5 st           | 7.5 lb      | 16"        | 17"        | 15" - 22"            | 15" - 22"        |
| UK-BTS118005 | Shower Chair with Arms - Aluminium     | 20"         | 31.5 st           | 6.8 lb      | 16"        | 17"        | 15" - 22"            | 15" - 22"        |