

Aspire Vida Folding Wheelchair - Self Propelled



- The Aspire Vida is a super lightweight transport wheelchair ideal for quick trips to the shops, appointments and catching up with friends.
- Ergonomic backrest and seat shape follow the contours of the user to maximise comfort and postural support.
- The 21 st safe working load caters to a broad range of users.
- Lightweight frame design ensures easy lifting in and out of cars, ideal for carers.
- Contoured and padded seating improves comfort and postural support.
- Attendant operated hand brakes with parking lock.
- Simple and intuitive folding for easy transportation.
- Removable leg rests and folding back for compact storage.
- Rear storage pocket for essentials.
- Available as attendant propelled or self-propelled.
- Anterior Seat Height - 19"
- Posterior Seat Height - 18"
- Anti-tippers available (WCA848275)
- Lap belt available (WCA848300)

Warranty Information

- 5 Years - Frame 1 Year - Fabric & Components

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Frame Weight	Backrest Height
UK-MWS449860	Aspire Vida Folding Wheelchair -	3'4"	27"	37"	21 st	25.8 lb	Orange	16"	18"	23.6 lb	17"

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Frame Weight	Backrest Height
	SP - Orange										
UK-MWS449865	Aspire Vida Folding Wheelchair - SP - Green	3'4"	27"	37"	21 st	25.8 lb	Green	16"	18"	23.6 lb	17"
UK-MWS449876	Aspire Vida Folding Wheelchair - SP - Mustard - 16"	41"	25"	37.5"	21 st	16.5 lb	Mustard	16.5"	16"	16.5 lb	17"
UK-MWS449888	Aspire Vida Folding Wheelchair - SP - Olive - 20"	41"	28"	37.5"	21 st	17.6 lb	Olive	16.5"	20"	16.5 lb	17"