

Aspire Bariatric Shower Commode



- Super Heavy Duty Bariatric design with 31" seat and 55 st SWL
- Anti-bacterial frame coating for improved infection control
- Removable arms to assist with positioning and allow lateral slide transfers
- Height adjustable, swing away footrests for optimal positioning and ease of transfers
- Padded seat and arms for maximum comfort
- 4x lockable castors maximise safety by providing carer access to a brake at all angles around the commode
- Comes standard with pan and carrier

Warranty Information

- Corrosion: 5 Years
- Components: 1 Year



Code	Display Name	Safe Working Load	Seat Width	Backrest Height
UK-BTC061200	Bariatric Shower Commode	55 st	31"	17"