

## Aspire Socialite Folding Wheelchair - Self Propelled



- The Aspire Socialite is an ultra light weight transport chair perfect for short trips and outings with friends.
- The generous 18 st safe working load caters to a broad range of users.
- Padded and breathable seat and back provide optimum comfort and postural support.
- Cut away arms improve table access and the offset uprights create hip space for coats and warmer clothing.
- 22" rear wheels for user self-propulsion.
- Sleek European design.
- Stylish contemporary colours.
- One of the lightest transport chairs available.
- Perfect for outings and appointments.
- Self operated hand brakes with parking lock.
- Simple and intuitive folding for easy transportation.
- Removable leg rests and folding back for compact storage.
- Contoured seating improves comfort and postural support.
- Rear storage pocket for essentials.
- Folded Dimensions (DxWxH) 26x15x29".
- Seat height anterior: 19 - 20".
- Seat height posterior: 17".

### Warranty Information

- 5 Years - Frame 1 Year - Fabric & Components

Code	Display Name	Total Length	Total Width	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Frame Weight	Backrest Height
UK-MWS449820	Aspire Socialite Folding Wheelchair - Self Propelled - Silver	41"	26.6"	18 st	28.7 lb	Silver	17"	17.7"	18 lb	16.1"
UK-MWS449830	Aspire Socialite Folding Wheelchair - Self Propelled - Red	41"	26.6"	18 st	28.7 lb	Red	17"	17.7"	18 lb	16.1"