

Aspire Evoke 2 HD Wheelchair



- The Evoke 2 HD is designed for the larger user who wants a lightweight, highly adjustable manual wheelchair.
- The Contoured Padded Overlay (shown in image #2) is a stylish recommended accessory for the Aspire Evoke 2. The contoured padded overlay offers the user increased comfort and is ideal for those who are in their wheelchair for longer periods of time. The mild contouring assists in positioning and pressure distribution, while the fabric allows air movement, increasing breathability. It is easily removed to allow washing when required
- Swing-in / out, height adjustable and removable foot rests – this facilitates user transfers and allows optimal positioning for the changing needs
- Strong, Nylon upholstery with a comfortable cushion and padded backrest. For increased comfort, a padded overlay is available as an accessory.
- Depth adjustable seating with rear angle adjustability – for more accurate positioning, customisability and increased user comfort.
- Quick release 24" solid rear tyres - the wheels are easily removed to reduce the weight of the chair for transport and storage.
- Crash Tested – to AS/NZ3616.196

Warranty Information

- 2 Years Frame
- 1 Year Components

Code	Display Name	Safe Working Load	Material	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
UK-MWS449540	Aspire Evoke 2 HD - 23" Wide	28 st	Aluminium	35.3 lb	Silver	16" - 22"	22"	17" - 19"	17" - 19"	17" - 19"
UK-MWS449550	Aspire Evoke 2 HD - 24" Wide	28 st	Aluminium	36.2 lb	Silver	16" - 22"	24"	17" - 19"	17" - 19"	17" - 19"