

Aspire OREGON Lift Recline Chair



- Simpler easy to use hand control function assisting cognition with only two buttons
- Contemporary medium size.
- Gentle tapered armrests suit the modern home environment.
- Integrated lift and recline functions for easy sit-to-stand.
- Generous padding in backrest for improved lumbar support.
- Headrest and armrest protectors included as standard.

Warranty Information

- 2 Years



Code	Display Name	Total Width	Total Height	Safe Working Load	Material	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
UK-CHP227150	Aspire OREGON Lift Recline Chair - Charcoal Fabric	31"	3'7"	23.5 st	Fabric	Charcoal	19"	19"	19"	19"	24"
UK-CHP227155	Aspire OREGON Lift Recline Chair - Latte Vinyl	31"	3'7"	23.5 st	Vinyl	Latte	19"	19"	19"	19"	24"