

## Aspire Rehab RX Advanced Tilt-In-Space Wheelchair



- ***The Rehab RX** – is designed for **the more complex user**, who requires assistance with **positioning, re-positioning and skin management**.*
- ***Tilt (26°) and backrest recline (30°)** – this facilitates weight re-distribution and the preservation of skin integrity. It also optimises user comfort, decreases user fatigue and assists with head control.*
- ***Contoured, built-up backrest and seat cushion** - the backrest has a gentle lumbar build up and moderate lateral supports; together with the abductor and adductor build-up at the seat cushion, midline positioning is optimised.*
- ***Advanced articulating legrests** – closely mimic the normal movement at the knee upon elevation, allowing the thigh to maintain contact with the seating surface, ensuring maximum support and skin contact with the cushion.*
- ***Breathability** - support surfaces offer multi-layered foam with breathable multi-stretch covers - for pressure re-distribution, comfort and to enhance skin microclimate.*
- ***Highly adjustable** - depth and width adjustable seat sizing. Height adjustable and removable armrests. Adjustable headrest.*
- ***Crash Tested** – to AS/NZ3616.196*

### Warranty Information

- 2 Years Frame
- 1 Year Components

Code	Display Name	Safe Working Load	Material	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
UK-MWS449710	Aspire Rehab RX - 17"-18" Wide	28 st	Steel	6.5 st	Silver	17" - 21"	16" - 18"	16" - 19"	16" - 19"	18" - 28"
UK-MWS449720	Aspire Rehab RX - 18"-20" Wide	28 st	Steel	7 st	Silver	17" - 21"	18" - 20"	16" - 19"	16" - 19"	18" - 28"
UK-MWS449740	Aspire Rehab RX - 20"-22" Wide	28 st	Steel	7.5 st	Silver	18" - 22"	20" - 22"	16" - 19"	16" - 19"	18" - 28"