

## Aspire Oslo Seat Walker



Sleek lines, elegant curves and eye-catching shades give the OLSO walker sophisticated poise. The aerospace inspired curved aluminium frame is highly durable and ultra lightweight for easy manoeuvrability, folding and lifting, while tool-free adjustments can suit individual user needs.

- **Ultra-Lightweight Design:** Constructed with aerospace-grade aluminium for exceptional strength, the 6.25 kg minimum weight makes it easy to lift and transport
- **Simple to Adjust:** Tool-free handle and backrest height adjustment accommodates users of varying stature, with push-button locking for secure customisation
- **Ergonomic Comfort:** Contoured, non-slip hand grips reduce wrist strain and the backrest is height and depth adjustable for greater support
- **All-Terrain Mobility:** 20 cm puncture-resistant PU wheels with swivel front casters and fixed rear wheels for smooth manoeuvrability indoors and outdoors
- **Portable & Compact:** One-click folding mechanism collapses the rollator into a slim profile (20 cm width) for effortless car boot storage or travel
- **User-Centric Accessories:** Secure zippered storage bag with 5kg capacity and built-in cane holder

## Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight
UK-WAF705320GR	Aspire Oslo Seat Walker - Green	26"	21"	33 - 37 "	21st	15 lb	18"	23"	23"	14 lb
UK-WAF705320PU	Aspire Oslo Seat Walker - Purple	26"	21"	33 - 37 "	21st	15 lb	18"	23"	23"	14 lb