

## Aspire Over Toilet Aid - Wide



- Safe Working Load (S.W.L) of 27.5 st. Tested to AS/NZS3973.7:2009
- Raised armrests – Assist the user transferring on and off the Over Toilet Aid
- Additional height increments on one leg – Accommodates uneven surfaces
- Height adjustable legs – Easy push button height adjustment for changing needs
- Comfortable, contoured seat with clip on attachment – For easy assembly and cleaning. Padded seat available, sold separately (BTT145110)
- All frames are individually treated for corrosion resistance
- Removable, moulded splash guard provided – Easily removed for transport and cleaning.
- Durable non-slip feet – Provide excellent grip and stability in wet conditions. Stainless steel washer within enhances durability
- Available in Aluminium or Treated Steel

### Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

Code	Display Name	Safe Working Load	Material	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
UK-BTT145010	Over Toilet Aid - Wide - Treated Steel	27.5 st	Treated Steel	9.3 lb	17"	21"	18" - 25"	18" - 25"
UK-BTT145015	Over Toilet Aid - Wide - Aluminium	27.5 st	Aluminium	7.5 lb	17"	21"	18" - 25"	18" - 25"