

Aspire Shower Stool



- Safe Working Load (SWL) of 27.5 st - Tested to AS/NZS39737:2009
- Clip on seat attachment - For easy assembly and cleaning
- Height adjustable legs - Easy push button height adjustment for changing needs
- Additional height increments on one leg - Accommodates uneven surfaces
- All frames are individually treated for corrosion resistance
- Raised armrests - Assist the user transferring on and off the chair
- Seat with anterior tilt - Facilitates standing from sitting and provides greater comfort for stiff or painful hips
- Raised armrests - Assist the user transferring on and off the stool
- Durable non-slip feet - Provide excellent grip and stability in wet conditions and a stainless steel washer within enhances durability
- Available in Aluminium or Steel

Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
UK-BTS119000	Shower Stool with Arms - Treated Steel	27.5 st	8.4 lb	12"	15"	18" - 26"	18" - 26"
UK-BTS119005	Shower Stool with Arms - Aluminium	27.5 st	6.8 lb	12"	15"	18" - 26"	18" - 26"