

Aspire Evoke 2 Wheelchair



- The Evoke 2 is designed for the user who wants a lightweight, highly adjustable manual wheelchair
- The [Contoured Padded Overlay](#) (shown in image #2) is a stylish recommended accessory for the Aspire Evoke 2. The contoured padded overlay offers the user increased comfort and is ideal for those who are in their wheelchair for longer periods of time. The mild contouring assists in positioning and pressure distribution, while the fabric allows air movement, increasing breathability. It is easily removed to allow washing when required
- Swing-in / out, height adjustable and removable foot rests – this facilitates user transfers and allows optimal positioning for the changing needs of the user
- Strong, Nylon upholstery with a comfortable cushion and padded backrest. For increased comfort, a padded overlay is available as an accessory
- Depth adjustable seating with rear angle adjustability – for more accurate positioning, customisability and increased user comfort
- Quick release 24" solid rear tyres - the wheels are easily removed to reduce the weight of the chair for transport and storage.
- Crash Tested – to AS/NZ3616.196

Prescribing Therapists: Click here to download our

[Clinical Assistance Guide](#)

Warranty Information

- 2 Years Frame
- 1 Year Components

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code: 05_122203104_0105_1_2

Wheelchair Manual –

Composite: 05_417_0105_1_2

Click [here](#) for more information

Compliance

AS3696-1:2008

