

Aspire Vogue Forearm Seat Walker / Rollator



- **Smart Posture System** – equipped with forearm supports that are height, depth and angle adjustable, the user can achieve optimal centre of gravity for improved balance and weight distribution.
- **Luxury look and feel** – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- **Lightweight and ergonomic design** with small turning circle for excellent manoeuvrability.
- **“AutoLocate”** - Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate.
- **Locks when folded** – makes transportation and storage easier.
- **Advanced locking handbrakes** – for increased user safety when walking and sitting.
- **Unbreakable TPE super soft wheels** – offers high durability and a smooth ride for increased comfort.
- *Unit weight includes accessories, bag and backstrap
- **Frame weight is minimal functional weight

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight
UK-WAF705400	Aspire Vogue Forearm	26"	24"	40" - 50"	23.5 st	1.5 st	White	18"	24"	24"	9.1 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Seat Width	Seat To Floor Height	Seat To Footrest	Frame Weight
------	--------------	--------------	-------------	--------------	-------------------	-------------	--------	------------	----------------------	------------------	--------------

Walker