

Freedom Bath Transfer Bench



- Removable backrest can be attached to either side of the bench
- Detachable armrest
- Height adjustable legs. Adjustment pins are spaced at 1' increments on three extension legs and 0.5' increments on the fourth to maximise stability
- Thicker, more durable powder coated aluminium frame to prevent scuffs and scratches
- Heavy duty sure grip suction rubber tips
- Purchased separately: Extension Legs (GSP650100)

Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
UK-BTS153100	Freedom Bath Transfer Bench - HBA424	20"	31"	32" - 38"	31.5 st	13.2 lb	16"	27"	19" - 25"	19" - 25"