

Aspire Shower Stool - Wide



- Reinforced wide frame – Provides stability for larger users
- Safe Working Load (S.W.L) of 27.5 st. Tested to AS/NZS3973.7:2009
- Contoured plastic seat with drainage holes and windows - Provides postural support, improves drainage and reduces skin adhesion
- Height adjustable legs – Easy push button height adjustment for changing needs
- Additional height increments on one leg – Accommodates uneven surfaces
- All frames are individually treated for corrosion resistance
- Raised armrests – Assist the user transferring on and off the elevated seat
- Durable non-slip feet – Provide excellent grip and stability in wet conditions and a stainless steel washer within enhances durability
- Available in Aluminium or Steel



Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

Code	Display Name	Safe Working Load	Material	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
UK-BTS119010	Shower Stool with Arms - Wide - Treated Steel	27.5 st	Steel	8.4 lb	14"	16"	17" - 24"	17" - 24"
UK-BTS119015	Shower Stool with Arms - Wide - Aluminium	27.5 st	Aluminium	6.8 lb	14"	16"	17" - 24"	17" - 24"