

Foot-up Drop Foot Orthosis



- The Foot-up is a lightweight ankle foot orthosis offering dynamic support for drop foot or similar complaints
- Provides visible improvement in gait by supporting the foot the moment it is raised
- Cushioned ankle wrap is made from a breathable, 3-layer material that can be worn comfortably for long periods of time
- Plastic inlay fits discreetly between the tongue and the laces of the shoe and is practically invisible when worn
- Sizing: Measure circumference just above the ankle bone

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information