

## Aspire Overbed Table - Recessed Tilt Top - Wide Base



- Internal width: 31.4" – can now be used around common Lift Recline Chairs
- Recessed Thermoform top to prevent spills
- Lightweight mobile tilting table positions overbed or chair
- Can be used for eating, drawing or other activities
- Angle adjustable to 45 degrees, ideal for homecare environment
- Side table does not tilt, suitable for drink and food
- Fully height adjustable
- "C" base shape fitted with four castors enables close positioning to the bedside or chair
- Stylish cream colour finish
- Powder coated mild steel frame construction

### Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Material
UK-BEA028100	Aspire Overbed Table - Recessed Tilt Top - Wide Base Thermoform - Cream	17"	31"	28" - 42"	Thermoform Top