

## Freedom Bath Transfer Bench with Rotating Seat- HBA425



- Soft, padded foldable backrest
- Soft, padded armrests
- Height adjustable legs with 0.5' increments on all legs
- Moulded seat with safety belt slides from end to end allowing easy patient transfer. Rotating seat allows easier access from any angle outside the tub
- 360° seat rotation lockable at 90o intervals
- Slip resistant suction rubber tips
- Lightweight aluminium construction

### Warranty Information

- 1 Year

freedom  
HEALTHCARE

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
UK-BTS153110	Freedom Bath Transfer Bench with Rotating Seat - 20 st	20"	31"	37" - 40"	21.5 st	19.4 lb	14"	17"	16" - 20"	16" - 20"