

Ceiling Hoist - Prism - Free Standing Gantry



prismmedicaluk
a prismhealthcare company

Features

- High-strength lightweight aluminium
- 200 kg max load
- 1.9 - 2.85 m adjustable track length
- 2.1 - 2.5 m adjustable height
- Stable 0.87 m
- Built-in easy-slide trolley
- Floor level adjusters
- Star Handle securing screws
- Safety latch to prevent uprights retracting
- Securing Bolts and lynch pins

Benefits

- Extends the time people can stay in their own homes
- No structural alterations required to the room or ceiling
- Complete safety, dignity and comfort for the client and carer
- Compatible with CP200 Portable Track Hoist
- Compatible with other portable hoists with QRS connectors
- Can be used to lift and transfer individuals from a bed, chair or similar fixture
- Easy to assemble, can be assembled by one person without tools (2 person assembly is recommended)

Warranty Information

- 1 Year

ndis

Item Code: 05_123612169_0105_1_2

Hoist – Composite: 05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

