

Aspire Overchair Table



- **Contoured Design:** Curved table top allows users to comfortably rest their arms whilst providing additional support to those with reduced upright sitting strength.
- **Personalised Width Adjustment:** Designed to accommodate recliners, chairs and beds up to 44" wide.
- **Raised Safety Trims** keeps items within arm's reach and prevents them from sliding off.
- **Height Adjustment** accommodates varying chair and bed heights.
- **Easy Clean Surfaces:** Durable Laminate table top and powder coated frame for fuss-free clean-up and usage.
- **Neutral Wood Grain:** Large table top accommodates food trays, personal items and provides a stable writing surface.
- **Locking Castors** for improved stability and effortless flexibility.

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour
UK-BEA028150	Aspire Overchair Table - Beech Laminate	36" - 44"	20"	28" - 36"	1.5 st	Beech