

Aspire Shower Chair - Wide



- Reinforced wide frame – Provides stability for larger users
- Safe Working Load (S.W.L) of 31 st. Tested to AS/NZS3973.7:2009
- Contoured plastic seat with drainage holes and windows - Provides postural support, improves drainage and reduces skin adhesion
- Height adjustable legs – Easy push button height adjustment for changing needs
- Additional height increments on one leg – Accommodates uneven surfaces
- All frames are individually treated for corrosion resistance
- Raised armrests – Assist the user transferring on and off the elevated seat
- Durable non-slip feet – Provide excellent grip and stability in wet conditions and a stainless steel washer within enhances durability
- Available in Aluminium or Steel

Warranty Information



- Corrosion: 5 Years
- Components: 1 Year

Code	Display Name	Total Width	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
UK-BTS118015	Shower Chair with Arms - Wide - Aluminium	22"	31.5 st	7.3 lb	16"	17"	15" - 22"	15" - 22"
UK-BTS118010	Shower Chair with Arms - Wide - Treated Steel	22"	31.5 st	8.6 lb	16"	17"	15" - 22"	15" - 22"