

Aspire Vogue Indoor Walker



- **Luxury look and feel** – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- Now supplied with bonus cane holder (ref. image)
- **Lightweight and ergonomic design** with small turning circle for excellent manoeuvrability.
- **No lifting & spillage** – meals can be collected and brought to the desired location without the need for lifting or carrying. The tray is also removable.
- **Sitting down for a meal** – when the user is sitting, the walker can be turned around and positioned comfortably in front of them, to enjoy their meal.
- **Unique folding design** – the slimline design reduces the footprint when folded but allows the folded walker to stand upright on its own.
- **Brakes** – the ergonomic ‘Soft Touch’ lockable design allows the walker to be secured in place with minimal hand pressure.
- **Removable storage bag** – utensils and other items can be stowed securely away, ideal for transporting personal belongings around the house.
- *Unit weight includes accessories, bag and tray
- **Frame weight is minimal functional weight

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour	Frame Weight
UK-WAF705100RE	Aspire Vogue Indoor Walker - Red/Black	25"	21"	34" - 38"	11.5 st	11.2 lb	Red/Black	4.4 kg

