



Vogue Carbon Seat Walker

Carbon Fibre Frame Extremely lightweight to maximise portability

One of the lightest seat walkers on the market. Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super-lightweight carbon fibre frame for a superior walking experience.

Less Than
5kg
11lbs



Deep Seat and Wide Backrest
For improved comfort and support



Compact Turning Circle
For improved mobility in tight spaces




1 YEAR WARRANTY



Sapphire
Short



Black
Medium



Emerald Green
Tall



Australia
1300 133 120 | aspirecare.com.au

United Kingdom
01622 541235 | aidacare.co.uk



Ergonomic Handles

These offer a comfortable grip during extended use. The flat-top design provides a stable resting place for arms while seated



Soft-touch Brakes

Lock the rear wheels to enhance safety and stability during use



Adjustable Handles

Allows for personalised positioning, helping to maintain proper posture



Simple Folding Mechanism

Locks for easy transportation and storage



Removable Storage

Zip bag secures up to 5kg of essentials



Accessories

With a cane holder included as standard, choose from our range of optional accessories

Size - Colour Code	Overall Length mm / "	Overall Width mm / "	Overall Height mm / "	Seat Width mm / "	Seat Height mm / "	Folded Width mm / "	Wheel Size mm / "	Minimum Functional Weight kg / lbs	Typical Functional Weight* kg / lbs
Short - Sapphire WAF705340	655	600	715 - 830	460	500	230	203	4.7	5.9
	26	24	28 - 32	18	20	9	8	10.4	13
Medium - Black WAF705350	655	615	725 - 940	460	550	230	203	4.9	6
	26	24	28 - 37	18	22	9	8	11	13.2
Tall - Emerald Green WAF705370	655	610	780 - 1030	460	620	230	203	5.2	6.3
	26	24	30 - 40	18	24	9	8	11.5	14

AS/NZ Approved | Vogue is registered in the ARTG, AUST R ARTG 271255 | *Includes accessories, bag, and backstrap.

Australia
1300 133 120 | aspirecare.com.au

United Kingdom
01622 541235 | aidacare.co.uk