



*Shown with pedals, not included as standard

Designed for the most demanding situations,
a solution for each department



1300 133 120 | aidacare.com.au

A MULTI-FUNCTIONAL WORK PLATFORM

The DELTA4 Ward Bed is the new generation of beds developed according to the needs of healthcare professionals. The DELTA4 bed becomes an active platform that simplifies and maximises the performances of healthcare professionals.

FEATURES

Bed exit alarm: Visual and acoustic alerts make the caregiver aware when a patient has left the bed (available on BEB050150 and BEB050155)

Predictive exit function: Alerts the caregiver that a patient may be attempting to leave the bed (available on BEB050150 and BEB050155)

Restless patient function: Signals the patient is experiencing restless behaviour so the carer can check and assist (available on BEB050150 and BEB050155)

Weight range function: Alerts carer that weight of patient which was preset has decreased, meaning they may have or may be attempting to leave the bed (available on BEB050150 and BEB050155)

INCENTRO System: the combined back and pelvic transition, reduces shear and friction forces, prevents slipping and keeps the body at a 30° angle, optimal for ventilation

Safe Exit Lights: 2 x LED lights signal whether the bed is at a safe exit height, promoting safe egress for the patient (orange is unsafe and green is safe)

Safe Exit control moves the platform to the minimum height and lifts the backrest up to facilitate the patient to get up from the bed

Auto Regression: backrest automatically shifts backwards thus increasing the space on the pelvic area by creating an extra space of 13cm

Biocote: metal components are epoxy powder coated and finished with an antibacterial additive for improved infection control

DeltaDrive: electric centre 5th wheel with integrated controls in headboard provides increased manoeuvrability for hospital staff (available on BEB050150 and BEB050145)

Weigh Scale: integrated system allows for precision weight measurements that can be set to automatically weigh a patient every 24 hours(available on BEB050150 and BEB050155)

4 sector safety sides with integrated electric controls which allows full autonomy of movement to the patient

CPR control on both sides allows the bed platform to be brought to the horizontal position immediately accessible

Mattress platform can be easily removed, all sections are designed for easy cleaning

Automatic self-brake once connected to the power net a sound alarm alerts the nursing staff about the unbraked bed and after 20 seconds the brake is automatically activated (available upon request)

X-ray tray together with HPL backrest allows extremely easy x-ray operations to be carried out without moving the patient (available upon request)



1 X-Ray Tray
(available upon request)

2 Side rails with integrated
electric controls

3 Safe exit

4 Removable mattress platform
for maintenance

5 High containment side
rails (42 cm from the
mattress platform)

6 Easily removable headboard and
footboard in case of emergency

7 Electric footrest

8 Integrated Weigh Scale
(available on BEB050150
and BEB050155)

9 Integrated 300mm
Bed Extension

10 Blanket and
'supervisor' holder

11 Led lights

12 Automatic Self Brake
(available upon request)

13 Electric Foot Pedals
(available upon request)

PRODUCT INFORMATION

- BEB050140** Delta 4 Ward Bed
- BEB050155** Delta 4 Ward Bed with Weigh Scale
- BEB050145** Delta 4 Ward Bed with Delta Drive
- BEB050150** Delta 4 Ward Bed with Weigh scale and Delta Drive

DELTA 4 beds come in four variants and have a range of accessories available to further improve their performance.



FEATURE	BEB050140	BEB050155	BEB050145	BEB050150
	Delta 4 Ward Bed	Delta 4 Ward Bed with Weigh Scale	Delta 4 Ward Bed with Delta Drive	Delta 4 Ward Bed with Weigh Scale & Delta Drive
Incentro System	✓	✓	✓	✓
150mm Castors	✓	✓	✓	✓
DeltaDrive Electric 5th Wheel	-	-	✓	✓
Integrated Weigh Scale	-	✓	-	✓
Bed Extension	✓	✓	✓	✓
Blanket Holder	✓	✓	✓	✓
X-Ray Tray (available on request)	-	-	-	-
4 Section Side Rails with Electric Controls	✓	✓	✓	✓
Underbed Light (available upon request)	-	-	-	-
Safe Exit	✓	✓	✓	✓
Removable Bed Ends	✓	✓	✓	✓
Removable Mattress Platform	✓	✓	✓	✓
CPR	✓	✓	✓	✓
Brake Alarm (Yellow/Green lights)	-	✓	✓	✓
Backlit Nurse Control Panel	✓	✓	✓	✓
Oxygen Bottle Holder (available upon request)	-	-	-	-
Radio Translucent Mattress Platform	✓	✓	✓	✓
Urine Bottle Holder	✓	✓	✓	✓
Accessory Bar Holder	✓	✓	✓	✓
Self Help Pole (BEA825090)	-	✓	✓	✓
IV Pole (BEA825101)	-	✓	✓	✓
Folding Monitor Tray (BEA825094)	✓	✓	✓	✓



1300 133 120 | aidacare.com.au



Exclusively distributed by