

MOBILITY & COMFORT



SCOOTERS

aspire HS268 Scooter UK-SCT671135



- Manual folding via hand lever
- Led headlights and brakelights
- Intuitive user control interface
- Approved for air travel
- Electro mechanical brakes
- Fast charging lithium ion battery

Length (mm)	Width x Height (mm)	Seat Width (mm)	Max Speed (mph)	Range (miles)	Product Weight (kg / st)	SWL (kg / st)
950 / 775 Folded	500 x 960 / 495 x 510 Folded	420	4	12.4	30.6 / 4.8 (incl. batteries)	115 / 18.1

aspire HS589 Scooter UK-SCT671270



- Led headlights and brakelights
- Front and rear wheel suspension
- Rearview mirror
- Adjustable swing away armrests
- Electro mechanical brakes

Length (mm)	Width x Height (mm)	Seat Width (mm)	Max Speed (mph)	Range (miles)	Product Weight (kg / st)	SWL (kg / st)
1290	625 x 1150	455	8	24.9	98 / 15.4 (incl. batteries)	136 / 21.4

aspire HS890 Scooter UK-SCT671288



- Electro mechanical brakes
- Swivel captain seat with seat sliding mechanism
- Hi-lo speed switch
- Rearview mirror & full lighting

Length (mm)	Width x Height (mm)	Seat Width (mm)	Max Speed (mph)	Range (miles)	Product Weight (kg / st)	SWL (kg / st)
1530	685 x 1245	480	8	23.3	114 / 18 (incl. batteries)	160 / 25.2

aspire HS898 Scooter UK-SCT671290



- Heavy duty hill climber
- Stylish mag wheels
- Led headlights, brakelights and indicators
- Front and rear wheel suspension
- Rearview mirrors
- Adjustable swing away armrests
- Electro mechanical brakes

Length (mm)	Width x Height (mm)	Seat Width (mm)	Max Speed (mph)	Range (miles)	Product Weight (kg / st)	SWL (kg / st)
1440	720 x 1360	520	8	24.9	143 / 22.5 (incl. batteries)	181 / 28.5

For full range of Mobility Scooters



LIFT RECLINE CHAIRS

aspire Aspire Oregon Single Action



- Tapered armrests suit the modern home
- Simultaneous adjustable back and leg support functions
- Lift and recline functions for easy sit-to-stand
- Generous padding in backrest for improved lumbar support

UK-CHP227150 Charcoal Fabric				UK-CHP227155 Latte Vinyl			
Size	Seat Height (mm)	Seat Depth (mm)	Seat Width (mm)	Backrest Height (mm)	Overall Height (mm)	Overall Width (mm)	SWL (kg / st)
Medium	480	480	480	620	1100	800	150 / 23.6

aspire Aspire Idaho Dual Action



- Tapered armrests suit the modern home
- Independent adjustable back and leg support functions
- Lift and recline functions for easy sit-to-stand
- Generous padding in backrest for improved lumbar support

UK-CHP227160 Charcoal Fabric				UK-CHP227165 Latte Vinyl			
Size	Seat Height (mm)	Seat Depth (mm)	Seat Width (mm)	Backrest Height (mm)	Overall Height (mm)	Overall Width (mm)	SWL (kg / st)
Medium	480	480	480	620	1100	800	150 / 23.6

aspire Aspire Raphael Quattro



- Raphael chairs are available in four sizes ensuring the perfect fit for all users
- Independent adjustable headrest, back, leg, and lumbar support functions
- Integrated footrest extension suits most user heights
- Integrated lift and recline functions
- Backlit hand control with USB charging point
- Adjustable headrest and armrest covers provided with microsuede and woven fabrics
- Comes with integrated all-day lithium battery backup

UK-CHP227301 Marina Grey				UK-CHP227300 Bronze			
Size	Seat Height (mm)	Seat Depth (mm)	Seat Width (mm)	Backrest Height (mm)	Overall Height (mm)	Overall Width (mm)	SWL (kg / st)
Medium	470	500	520	580	1090	780	130 / 20.5

For full range of Electric Lift Recline Chairs



BEDROOM

 **aspire**
ComfiMotion® Collection



Aspire Comfimation Activ Care



- Fully adjustable hi-lo bed with Trendelenburg
- Ultra low height mattress platform of 195mm
- USB charger for devices
- Night assist light for safe entry and exit
- Box-on-box design facilitates easier entry and exit with true flush alignment
- Personalised positioning memory
- Press and hold safety handset
- Massage function for neck, shoulders and back

Size	Code	Overall Length (mm)	Overall Width (mm)	Overall Height (mm)	SWL (kg / st)
Single	UK-BEB046740	2030	920	195 – 700	250 / 39.4
King Single	UK-BEB046750	2030	1070	195 – 700	250 / 39.4

AIDACARE LIMITED

26 Kings Hill Ave
Kings Hill
West Malling
Kent
ME19 4AE

01622 541235 | aidacare.co.uk

For full Comfimation range visit

