




OSLO SEAT WALKER

Sleek lines, elegant curves and eye-catching shades give the OSLO seat walker sophisticated poise. The aerospace grade aluminium frame is highly durable and ultra lightweight for easy manoeuvrability, folding and lifting, while tool-free adjustments can suit individual user needs.



- **Ultra-lightweight Design:** Constructed with aerospace-grade aluminium for exceptional strength, the 6.25kg minimum weight makes it easy to lift and transport
- **Simple to Adjust:** Tool-free handle and backrest height adjustment accommodates users of varying stature, with push-button locking for secure customisation
- **Ergonomic Comfort:** Contoured, non-slip hand grips to reduce wrist strain and adjustable backrest for tailored support
- **All-terrain Mobility:** Puncture-resistant PU wheels with swivel front castors and fixed rear wheels improve manoeuvrability
- **Portable & Compact:** One-click folding mechanism collapses the rollator into a slim profile (20cm width) for effortless car boot storage or travel
- **User-centric Accessories:** Secure zippered storage bag with 5kg capacity and built-in cane holder



Fold
Pull handle up to collapse & lock



Brake
Lock: down
Unlock: up



Adjust
Press & hold to adjust

SAPPHIRE BLUE

EMERALD GREEN

1 YEAR WARRANTY
Components

5 YEAR WARRANTY
Aluminium Frame



| Code | Description | Overall Height | Overall Length | Overall Width | Seat Depth | Seat Width | Seat Height | Folded Width | Wheel Size Front & Rear | Min Functional Weight | Typical Functional Weight |
|-------------|---------------|----------------|----------------|---------------|------------|------------|-------------|--------------|-------------------------|-----------------------|---------------------------|
| WAF705320BL | Sapphire Blue | 33" - 37" | 26" | 21" | 8" | 18" | 23" | 8" | 9" & 8" | 6.3kg | 6.8kg |
| WAF705320GR | Emerald Green | 33" - 37" | 26" | 21" | 8" | 18" | 23" | 8" | 9" & 8" | 6.3kg | 6.8kg |